

OPTIMIZATION

EFFICIENCY



LIQUID SCINTILLATION COUNTING : SMEARS AND AQUEOUS SAMPLES FROM INDUSTRIES

LSC MEASUREMENTS

OPTIMIZATION OF PARAMETERS

CONSUMABLES, SAMPLING, METHODS

CHOICE OF CONSUMABLES

Liquid Scintillation Counting (LSC) is the detection of a radiation through the scintillation light produced in solutions. It has the advantage of mixing the sample within the sensitive medium, allowing the detection of low-energy radiations.

The choice of consumables has a strong impact on the quality of countings made by liquid scintillation. However, optimization of measurement conditions consists in maximizing the Figure of Merit :

$$FoM = \varepsilon / b^2$$


Where ε is the detection efficiency and b the background noise level.

Based on this criteria, the most suitable consumables are 20 mL polyethylene vials, containing:

- 10 mL of UltimaGold scintillation liquid to measure radioactivity on smears with cellulose filter folded and rolled, placed at the bottom of the vial
- 10 mL of UltimaGold scintillation liquid to measure radioactivity in liquids with 2 mL of liquid sample mixed.

MEASUREMENTS PARAMETERS

Radionuclides that can be searched by LSC are for example ^3H , ^{14}C , ^{137}Cs or ^{40}K . Their energy spectra being different, definition of counting windows is also tuned on the basis of the FoM maximization.



Under these requirements, best counting windows suited to the desire radionuclides are :

- [0:8] keV for ^3H ;
- [12:80] keV for ^{14}C (picture above);
- [25:400] keV for ^{137}Cs ;
- [100:750] keV for ^{40}K .

EFFICIENCY AND DETECTION LIMITS

The detection limit is ~100 Bq/L for ^3H in a liquid sample and ~50 Bq/L for ^{14}C , ^{137}Cs or ^{40}K considering ten minutes of counting time.

This difference is mainly due to the detection efficiencies associated to these radionuclides which are ~40% for ^3H and ~75% for ^{14}C , ^{137}Cs or ^{40}K and background which is prevalent at low energies.

Find more in "Optimization of liquid scintillation measurements applied to smears and aqueous samples collected in industrial environments", A. Chapon et al., Results in Physics 6 (2016) 50–58.



Arnaud Chapon
ATRON Scientific and technical manager

ATRON - 14 allée des Vindits
50130 Cherbourg-en Cotentin, France
+33(0)2.61.81.99.90
achapon@atron.fr



Gilbert Pigrée
Radiation protection engineer

IMOGÈRE - University of Caen, Normandie
B.P. 5186, 14032 Caen cedex
+33(0)2.31.56.63.90
gibert.pigrée@unicaen.fr